

# Frontiers in Science

## Public Lecture Series

Sponsored by the Fellows of Los Alamos National Laboratory

### Measuring Small Magnetic Fields in Living Systems: From Understanding the Brain to Cancer Detection


Dr. Michelle A. Espy, Biophysics Group

Española: Wednesday, 2 October 2002, at 7:30 PM  
Center for the Arts Theater, Northern New Mexico Community College

Santa Fe: Wednesday, 9 October 2002, at 7:30 PM  
James A. Little Theater, New Mexico School for the Deaf

Los Alamos: Thursday, 10 October 2002, at 7:30 PM  
Duane W. Smith Auditorium, Los Alamos High School

Taos: Wednesday, 16 October 2002, at 7:30 PM  
Taos Convention Center



### The Complexity, Simplicity, and Unity of Living Systems

Dr. Geoffrey B. West, Elementary Particle and Field Theory Group

Santa Fe: Wednesday, 13 November 2002, at 7:30 PM  
James A. Little Theater, New Mexico School for the Deaf

Taos: Thursday, 14 November 2002, at 7:30 PM  
Taos Convention Center

Española: Wednesday, 20 November 2002, at 7:30 PM  
Center for the Arts Theater,  
Northern New Mexico Community College

Los Alamos: Thursday, 21 November 2002, at 7:30 PM  
Duane W. Smith Auditorium, Los Alamos High School



Admission Is Free

<http://stb.lanl.gov/fellows/fellows.html>

Los Alamos  
NATIONAL LABORATORY


# Measuring Small Magnetic Fields in Living Systems: from Understanding the Brain to Cancer Detection.

---

**Michelle Espy**  
**SQUID Team**


**Biophysics Group, P-21**

**Los Alamos National Laboratory**



**The Superconducting Quantum Interference Device (SQUID) is the world's most sensitive detector of magnetic fields.**


**This talk is about how SQUIDs are used to measure the tiny magnetic fields associated with the processes of life in health and illness.**



# Your body is electric, and thus magnetic



Neuron in the human brain



Human heart

Your body and brain are full of nerves which transmit and receive information electrically, producing magnetic fields that provide a window into their function.


## Even if the process we are interested in isn't magnetic, sometimes we can make it magnetic

---

Take tiny magnetic beads, smaller than a blood cell, and use them to label cancer.

This provides you with a tool for looking for cancer by looking for where the magnetic signal is.

You can even try to use these magnetic beads to kill the cancer.



**Magnetic particle bound to cell**



# The SQUID team

---

Robert Kraus, Michelle Espy, Andrei Matlachov, Val Armijo, Shaun Newman, Dave Clark, Carl Kumaradas, Petr Volegov, Mark Peters

Students: Anne Chamberlin (LAHS/Lewis and Clark), Shawn Forrest (U. Utah), Jessica Lamb (Dartmouth/Cornell), Jonatan Mattson (U. Chalmers, Sweden), Walter Roybal (UNM/Pojoaque). We are ALWAYS looking for students!


External Collaborations: UNM, VA Hospital, U. Nebraska, MIND Inst. etc.



Biggest magnetic  
fields (NHMFL)

# Magnetic Fields

This scale is logarithmic  
(like the Richter scale  
for earthquakes) so each  
division indicates a  
magnetic field 10 times  
larger than the one  
below it!



The sensitivity is both  
the best and worst things  
about a SQUID – and  
why they aren't very  
common

# If SQUIDs are so great, why aren't they everywhere?

---

Some systems mentioned tonight have electrical counterparts (EKG for hearts, EEG for brains) that are more common. Why?



Because SQUID systems are harder to build.





# Superconducting Quantum Interference Device

---



The devices are ~1" in size – but most of that is packaging – places to connect to.

Actual SQUID sizes are on the order of 10-100 microns (like a human hair).

They are made very much like regular computer chips, but with superconducting material.

# SQUIDs are tools

---



The SQUID team is an instrumentation team: we build systems that use SQUIDs as tools to measure the magnetic fields associated with processes of interest:


The heart – function when healthy, or diseased.



The human brain – localizing function for surgery, epilepsy, effects of illness, how the brain is organized and processes information.



Cancer – early detection and diagnosis, better treatments.



Bioassay- Looking for DNA sequences, viruses, pathogens.

# Measuring Small Magnetic Fields in Living Systems: from Understanding the Brain to Cancer Detection.

---

Magnetic fields provide the clues to the biological processes.

SQUIDs can measure these fields totally non-invasively:


nothing touches you

no radiation


not loud

doesn't hurt


no injections



# Magnetic Fields: they are all around us




Bar magnet (refrigerator magnet),  
Strong on surface, but goes away  
quickly with distance.



Magnetic field of the earth.  
10,000 smaller than a  
refrigerator magnet, but we see  
it's effects: i.e. the compass.




MRI systems,  $\sim 10$  times stronger than  
refrigerator magnet, over a large area.  
Produced by large currents. Often the  
magnets are superconducting, meaning that  
they are kept VERY cold and the big currents  
flow without resistance.




Current will produce a  
magnetic field that is  
proportional to its size.  
CD player power adapter cord  
at 10 inches has field a million  
times smaller than the field of  
refrigerator magnet.


National High Magnetic Field  
Laboratory (Los Alamos, Florida  
State, U. Florida) has magnets  
that produce fields 100 – 1000  
times stronger than the fields  
from a refrigerator magnet





# Small Magnetic Fields in Living Systems: really small!



Magnetic field from electrical activity in the heart (a billion times smaller than the refrigerator magnet's field)



Magnetic field from electrical activity in the brain (roughly a trillion times smaller than the magnetic field from your refrigerator magnet)



Tiny magnetic particles can be attached to an agent that seeks a certain type of cell receptor, i.e. cancer. This gives the cancer a magnetic label about a billion times smaller than the magnetic field from a refrigerator magnet when you measure on the surface of the body.





# Measurements with SQUIDs

---

We use small passive detectors called SQUIDs, superconducting quantum interference devices, to detect these magnetic fields.


The SQUID is the most sensitive detector of magnetic fields known, capable of measuring magnetic fields billions of times smaller than the magnetic field of the earth.






# Outline

---





- Introduction to the Biophysics Group
- Introduction to SQUIDs
- What we do with SQUIDs (biological applications)
  - Heart function
  - Brain Function
  - Magnetic Cancer detection and treatment
  - Bioassay


# The Biophysics Group




Protein Crystallography




SQUID applications




Biomolecule detection



Multi-modality brain imaging



Computational Neuroscience




High Speed Electronics

- **Who we are**

- ***Physicists, Physiologists, Neuroscientists, Mathematicians, Technicians, Computer Scientists, Electrical Engineers, Support Staff, Students....***

- **Our Mission**

- ***Contribute to the understanding of biological phenomena by means of the scientific, technical, and conceptual resources of Physics***



---

**What is a SQUID anyway?**

# S is for Superconducting


A SQUID is made up of a loop of superconductor: current (electrons) flow with no energy loss or resistance.



In normal metal electrons bump into atoms and lose energy as they flow (resistance – or why a light bulb lights). In a superconductor the electrons don't bump into anything, and current flows with no loss.

But a SQUID is a superconductor with a twist (or two).

A SQUID has two spots that are not good superconductors. They are like kinks in a hose – so at some current the electrons can't get through without bumping. And you get resistance, which can be measured as a voltage across the SQUID.





# **Q is for Quantum I is for Interference**


---

**The electrons are also like waves (that's Quantum Mechanics). So you have the same effect from the two “kinks” in the SQUID that you have dropping two pebbles into water – wave interference.**

**The waves also depend on the magnetic field through the SQUID. This makes the voltage change depending on the magnetic field. SQUIDS can detect changes in magnetic field a million times smaller than the magnetic field of the earth passing through an area the size of a red blood cell.**



## D is for device: the basics of a SQUID system



All SQUIDs must be cold, and are usually cooled with either liquid helium (- 452 F) or liquid nitrogen (- 321 F).


The temperature of the surface of Pluto is - 378 to - 396 F




# Functional Heart Imaging with Magnetocardiography (MCG)

---

The heart is electrical and the magnetic fields we measure have the same source as electrocardiography: ECG or EKG



Magnetic field map above a beating heart



QRS complex of a healthy heart

# Functional Heart Imaging with Magnetocardiography (MCG)

---


EKG is more common and electrodes are less susceptible to problems with external magnetic signals. Why do MCG instead of EKG?

Non-contact, less preparation (mass screening, burns)

Complementary information, no signal losses in tissue

May be more sensitive to certain conditions

Fetal MCG – there is no good electrical way to measure the QRS of the fetal heart






EFMT, Bochum



with permission: CTF, Canada

At Los Alamos we have used our computational resources to help solve problems like removing the mother's heart beat from a fetal MCG signal.

# MagnetoCardiography (MCG)




liquid nitrogen



liquid helium

Frog MCG



# Functional Brain Imaging with Magnetoencephalography (MEG)


---

Use magnetic field information to probe where function is located and the timing of function. Why is this type of imaging important?

**Study Brain Function:** How, When, Where does the brain respond to stimuli?

**Diagnostic Tool:** Diagnose epileptic spikes and locate where they are initiated. Locate a region of brain injury.

**Pre-surgical Localization:** 'Map' primary brain centers prior to surgery.



# Functional Brain Imaging with Magnetoencephalography (MEG)


MEG/EEG measure responses directly from neurons. MEG is similar to EEG, but less common. Which is better? Ideally you would do both.

There are other methods to study brain function.

Implanted electrodes (invasive).


PET (radioactive)  
and fMRI, like an MRI

PET and fMRI only tell you about blood flow, which is slower ( $\sim 1$  sec) than neuronal activity ( $1/1000$  sec).




# What good is the timing information?

- Sources of epilepsy





- The effects of medication on the brain
- How the brain recovers from injury
- How the brain is processing information






# Early ideas of how the brain worked....



Planning ahead, prediction, speech  
(lower mostly left), motor (thin strip in  
front of central sulcus)


Receives nerve impulses for  
pain, temperature, touch, and  
pressure: primary sensory area.





Visual information



Hearing, memory processing, may be  
involved in integration of multiple  
sensory functions i.e. hearing, and touch.

©EnchantedLearning.com







# The Brain is electrical, and thus magnetic




The brain is full of nerve cells



The nerve cell receives information, in the form of electrical signals, from other cells through connections or synapses on its dendrites and sends information by way of its axon.




Magnetic Field Line






# Magnetoencephalography MEG



SQUIDs measure magnetic fields outside of the head due to neuronal activity



# Mapping Brain Function... and more

---

- Where function is located
  - pre-surgical, site of epilepsy
- The time course of activity
  - faster than 1/1000 of a second responses directly from neurons
- To understand how the brain processes information and what this means to healthy brains and in the treatment of disease



LANL MEG System

Totally passive measurement, no noise, radiation, head restraints, etc.






# The Los Alamos MEG system

---



MEG systems are available commercially – but ours is different

Most systems use a series of complicated coils and expensive shielded rooms to reduce background fields.

Our system has a superconducting shield built in – which means you don't need the coils or as expensive of a shielded room



# The Los Alamos MEG system




Corian®-like support and mounting structure.  
(DuPont/LANL Patent)

The LANL Whole-Head SIS-MEG System

Lead Superconducting Imaging Surface (SIS) "helmet".


## Superconducting Image-Surface Concept




Left: Traditional Gradiometer coupled to SQUID.

Right: Superconducting Imaging-Surface Concept showing source & source image ( $M$  and  $M_{img}$ ).

# The Los Alamos MEG system



Corian®-like support and mounting structure.  
(DuPont/LANL Patent)



The LANL Whole-Head SIS-MEG System



Lead Superconducting Imaging Surface (SIS) "helmet".

The cryogenic column is immersed in liquid helium at 4K








# From raw data to results




Computation




# Calibration with precision phantoms




Experimental data




Model data after CS transf.  $R_z[16.9] T[-0.2 \ 0 \ -1.1]$



# From raw data to results



LANL SIS Sensors  
Somatosensory Localized Current Dipole N20






# Where To Next?

---



- **System is at the VA Hospital in Albuquerque for performance studies**
  - The VA hospital has a commercial system that they use for epilepsy, pre-surgery, and other studies.
  - We want to compare our system to theirs
  - first tests with both phantoms and people
  - Our system will be as good (better) and need less shielding
- **At present our spatial resolution of phantom sources is 1/100 of an inch (better than the best reports in the literature)**
- **We are extremely pleased!**



# Why are we doing MCG or MEG?


---

- **Systems exist commercially, but at Los Alamos we have the resources and to try and push the technology in new directions.**



# Bioassay and MagnetoCarcinotherapy: labelling what you are interested in with tiny magnetic particles

---

- In collaboration with
  - University of Nebraska, University of New Mexico, Bioscience Division, Spallation Neutron Source Division, Chemistry Division
- Small (bacteria sized) magnetic particles (labels) attached to agents that seek out and label specific biological targets
- SQUIDs to detect and/or image particles



**The material must be coated, sorted, and attached to the right agents.**




A magnet made of this material is 10 times stronger than your average refrigerator magnet of the same size.


**This step of attaching the labels to something that is specific only to a certain type of cell (cancer) or DNA fragment, etc. is an active area of research – its called molecular targeting.**

**Both of the next two projects I will discuss are VERY new and involve using these magnetic labels to monitor biological processes.**

# Carcinotherapy: Target the cancer




# Carcinotherapy: Image




- Imaging of magnetic particles at tumor can be done with same techniques used in MEG

Visualization



Modeling & Computation



Instrumentation




*Localizing concentrations of magnetic particles *in vivo* uses techniques developed for functional brain mapping*



- To image a tumor of 1/100 to 1/2 inch in size, 1 inch away, imagining you can deliver 20x more material to the tumor than to healthy tissue, you need ~20 nanograms (about 1/2000 of a grain of salt).



# Carcinotherapy: Treat




- The needle on a compass wants to line up with a magnetic field.
- Similarly, our magnetic particles are like tiny compass needles and will line up in a magnetic field.
- If the magnetic field rotates, our magnetic bead will follow it. Acting like a little stirring bar.



# Carcinotherapy: Treat

- The rotating magnetic heats particles  $>45$  F in  $\ll 1$  min
- Tumor is killed with minimal damage to nearby healthy tissue



Hardware  
supplied  
by SNS  
Division



Put a sample (tissue, salt water, etc.) inside this box and see how the heating progresses in the material with magnetic beads vs. the material without.

# Carcinotherapy: Treat

---

- To treat, one needs ~ 5 micrograms (about 1/10 a grain of salt) in a tumor of size ~1/100 inch.
- So far we are only testing on samples.
- But the idea is that if you were applying this treatment to a person, the tumor would be killed and when the treatment was over the body would eliminate the beads through the kidneys.
- The benefit to this method is that it targets primarily **ONLY** the cells that are unwanted.
- The idea of heating tissue with magnetic fields to kill cancer is decades old and called Hyperthermia
- Our idea is different because the magnetic material is so strong. But targeting is still being developed!




# Bioassay: What is it?

---

Bioassay is a way of determining the strength or biological activity of a substance. This can be done on a living system or on a sample (tissue, blood, etc.).

For example, lets say you wanted to know if a person had a certain DNA sequence, or a certain virus...

Desire to look for a large number of agents simultaneously



# Bioassay: Why?

---

**Imagine you are some sort of bio-detective:**

**What is in that envelope of white powder?**

**What strain of flu does the patient have?**

**Maybe you don't really know what you are looking for and want to match your evidence with many possibilities....**




# General Approach

**Attach a label to an agent specific for the target**

**Let the label+agent mix with the sample (agent will bind specifically with the target)- in vitro we do a 100 to 1000 times better than in the body!**


**Targets are made such that target + tag is fluorescent – shining a laser and looking for the light can tell if you've got a match**






# Look for more than one thing

Use different labels:



1) Colored beads with lasers sorting, Flow Cytometry  
The National Flow Cytometry Resource is at LANL

2) Or you could use magnetic beads (little bar magnets) with different strengths




# Bioassay with magnetic nanoparticles


---

## What is a magnetic nano-particle?

A magnetic bead that is **VERY** small, ranging in sizes from the width of human hair to smaller than a bacteria. The beads will have different magnetic strengths (sizes) but are coated so they are roughly the same size on the outside.



# Bioassay with magnetic nanoparticles: First you sort




# Bioassay with magnetic nanoparticles: Incubate


---

Attach a label to an agent specific for the target

Let the label+agent mix with the sample

Targets are made such that target + tag is fluorescent – shining a laser and looking for the light can tell if you've got a match






**Sample (in tube) introduced in flow**



**Illuminated by laser,  
fluorescence tells you there was  
binding, don't really do this  
with your eye!**

**Flows past SQUID to tell you what  
kind of label**

# Bioassay with magnetic nanoparticles: present status




Signal from a bead passing under two SQUIDs. By knowing the separation and the timing shift we can calculate that the bead is traveling about 20 inch/sec.






## Concluding remarks



- SQUIDs are among the world's most sensitive magnetic field detectors.
- The Biophysics group applies SQUIDs to problems in support of the laboratory's mission, from brain imaging to basic physics.
- Never use a SQUID unless you have to. (Weinstock's Law)
- When nothing else will work- try a SQUID!




# Concluding remarks

---

**The brain and heart are electric, and thus magnetic.**


**We use SQUIDs to measure these magnetic signals, to understand how the brain and heart function, when people are healthy or suffering from an illness.**



**One can also make the processes of interest magnetic and use the sensitivity of SQUIDs to monitor them. This area is new and we are just beginning.**

**At Los Alamos we are building several systems to measure magnetic fields more effectively.**

**SQUIDs are the tool, but the goal is to improve the quality of life.**





Thank you!